



PRODUCT SPECIFICATION



Product: SALVEQUICK TEXTILE (FABRIC) XL
ADHESIVE BANDAGE (PLASTER)
Reference code: 6470



1. General description

Adhesive bandages of two sizes, small and big, packed in a cardboard refill.
The refill fits into the Cederroth (Salvequick) dispenser for adhesive bandages.

The adhesive bandage consists of a backing, evenly spread with an adhesive coating. An absorbent pad is placed centrally on the dressing. The adhesive is covered with a release-coated foil that consists of two parts partly overlapping each other. When the adhesive bandage is pulled out from the package, the release-coated plastic foil remains behind. One half of the dressing can now be applied on the skin using one hand only. The rest of the coated paper can then be removed and the adhesive bandage is ready for use.

2. Intended use

The adhesive bandages are intended to be used as a mechanical barrier and for absorption of exudates from minor wounds.

3. Technical description

Product dimensions

Small	Adhesive bandage:	80 x 30 mm
	Wound pad:	24 x 30 mm
Big	Adhesive bandage:	80 x 60 mm
	Wound pad:	24 x 60 mm

Material

Backing:	Elastic, beige fabric made of viscose/polyamide with hydrophobic treatment
Adhesive:	Polyacrylate adhesive
Woundpad:	Nonwoven viscose/polypropylene/polyethylene covered with a polyethylene wound contact layer
Facing paper:	Single sided siliconised release coated paper/foil
Wrapping paper:	Plain coldseal paper



PRODUCT SPECIFICATION



4. Composition

Refill

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
45164700	Small textile adhesive bandage	14
45164701	Big textile adhesive bandage	7
Dimensions (l x w x h) mm:		30 x 87 x 100

Consumer package

Number of refills in consumer package:	6
Dimensions (l x w x h) mm:	100 x 150 x 90

Transport package

Number of consumer packages in transport package:	10
Dimensions (l x w x h) mm:	520 x 160 x 190
Weight (kg):	2,5

5. Labeling

Product (refill)

Cederroth contact information
Product name, Reference code, CE, Lot-number,
Expiry date, EAN-code, Do not reuse, Size,
Quantity

Consumer package

Cederroth contact information
Product name, Reference code, CE,
Lot-number, Expiry date, EAN-code
Quantity

Transport package

Cederroth contact information
Product name, Reference code, CE,
Lot-number, Expiry date, EAN-code, Quantity

6. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

7. Accessories

The following accessories can be used with the product but are not necessary.

<u>Ref #</u>	<u>Product name</u>
490700	Salvequick Plaster Dispenser



Authorization:

QA&RA Mgr
2015.08.20 14:21:43 +02'00'